

Access Troubleshooting Guide

General and Positioning

- Check if the device power is on
- If applicable, check if the other equipment or devices (such as switch, mouse, camera, joystick etc.) are properly connected
- Check if the access method type is properly enabled or selected within the software application being used
- **Positioning**
 - i** For a permanent and optimal solution, consult with an occupational and/or physical therapist for assistance with maximizing function in various body positions
 - i** Always remember to **position device to user, NOT position user to device!**

Problem	Solutions
Unable to make selections on device and maintain position of body	As a temporary solution , stabilize the trunk using pillows or towels to see if this improves functional use of device Try a different access method
Uses device in various positions (e.g. wheelchair, bed, chair, etc.)	Evaluate access methods and mounting for each of the individual's positions
Fatigues quickly in current position	Re-evaluate and/or adjust the current position Try a different access method

Eye Tracking

- For difficulty of any kind, check position of device and person. Reposition if necessary
- Use Track Status Box for calibration

Problem	Solutions
Eyes move better vertically	Reposition device up or down Use device with a smaller screen Adjust calibration area to the part that the user is strongest on
Eyes move better horizontally	Reposition device left or right Use device with a smaller screen Adjust calibration area to the part that the user is strongest on

Calibration results are poor	<p>Recalibrate and/or reposition</p> <p>Change calibration target (color, type, animation) and/or background color - this may increase individual's attention to the task</p> <p>Decrease number of calibration targets - this will help with decreased attention and fatigue</p>
Accuracy seems to get worse with continued use	<p>Look away or close eyes for a moment</p> <p>Gaze at border of target, just above or just below target</p> <p>Check track status to make sure camera is still picking up eyes</p>
Glare from glasses interferes with use	<p>Clean glasses</p> <p>Reposition device (raise and tilt device more toward eyes)</p> <p>Limit or eliminate light coming from behind device towards individual's eyes</p> <p>Consider dark frames instead of wire frames</p>
When using dwell to select, difficulty holding eye gaze on buttons for selections	<p>Shorten dwell time</p> <p>Change highlight to invert rather than outline</p> <p>Try other options such as switch or blink</p>
Selections happen too fast or accidentally	<p>Lengthen dwell time</p> <p>Use Pause feature</p> <p>Change highlight to outline rather than invert</p> <p>Try other options such as switch or blink</p> <p>Try gaze feedback options such as shrinking dot or clock</p> <p>Provide a rest area visual on the device or screen</p>

Mouse

Problem	Solutions
Difficulty moving the mouse to all areas of the screen	<p>Adjust mouse speed on device</p> <ul style="list-style-type: none"> • Faster speed requires more control but less range of motion • Slower speed requires more range of motion but less control) <p>Move mouse into different positions (e.g. on the lap, closer, farther, right left)</p> <p>** For head-controlled mouse:</p> <ul style="list-style-type: none"> • Check positioning of person and device • Teach how to re-center mouse in work area <p>Try a different mouse</p>

Difficulty using mouse button to click	<p>Change to dwell for selection (hold cursor over button to select versus clicking)</p> <p>Use external switch for selection</p> <p>Try a different mouse</p>
When using dwell for selection, unable to hold mouse still long enough for selection to occur	<p>Shorten dwell time</p> <p>Change the highlight to invert rather than outline</p> <p>Add auditory feedback</p>
When using dwell for selection, makes accidental selections	<p>Lengthen dwell time</p> <p>Move mouse in different position (e.g. on the lap, closer, farther, right, left)</p> <p>Use switch for selection</p>
(For head mouse) Cannot follow mouse pointer on screen	<p>Modify mouse pointer properties (color, size, shape) if the options are available</p> <p>Change highlight to invert rather than outline</p>

Scanning

Problem	Solutions
Cannot activate switch accurately, consistently or in a timely manner	<p>Check switch position</p> <p>Try a different switch</p> <p>Try a different body part (hand, finger, head, shoulder, elbow, knee, foot, isolated muscle movement)</p>
Hits switch too early	<p>Check switch placement</p> <p>Increase scanning speed</p> <p>Add auditory feedback</p> <p>Adjust hold time</p> <ul style="list-style-type: none"> • If user hits too early, increase hold time • If user cannot hold switch down, decrease hold time <p>Adjust release time</p> <ul style="list-style-type: none"> • If user hits switch multiple times in a row on accident, increase release time • If user cannot select when activating switch multiple times, decrease release times <p>Try different scan pattern or method</p>

Hits switch too late	<p>Check switch placement</p> <p>Decrease scanning speed</p> <p>Add auditory feedback</p> <p>Adjust hold time</p> <ul style="list-style-type: none"> If user cannot hold switch down, decrease hold time <p>Adjust release time</p> <ul style="list-style-type: none"> If user hits switch multiple times in a row on accident, increase release time If user cannot select when activating switch multiple times, decrease release times <p>Try different scan pattern or method</p>
Takes too long to scan and make selection	<p>Increase scanning speed</p> <p>Change scan pattern</p>
Difficulty making accurate selections	<p>Try the zoom / magnify option</p> <p>Try different scan pattern or method</p>
Loses attention to task	<p>Adjust scan pattern</p> <p>Change highlight to invert rather than outline</p> <p>Add auditory feedback</p>

Touch

Problem	Solutions
Difficulty making accurate selections	<p>Increased hold time</p> <p>Increase release time</p> <p>Change Touch selection type (Touch Enter, Touch Exit)</p> <p>Try a keyguard</p>
Cannot access entire screen	<p>Check device position</p> <ul style="list-style-type: none"> Closer requires more motion but less control Farther requires less motion but more control <p>Try another access method</p>